



ALCON/CIBA Contact lens parameters

Alcon/CibaVision Daily disposable contact lenses

Contact Lens Name	Contact Lens Specifications	Contact Lens Parameters
Dailies Total 1	Base Curve: 8.5 Diameter: 14.1 Dk/t @-3.00D 156	-0.50 to -6.00 (0.25D steps) -6.50 to -10.00 (0.50D steps)
Focus Dailies AquaComfort Plus	Base Curve: 8.7 Diameter: 14.0 Dk/t @-3.00D 26	+0.50 to +6.00 (0.25D steps) -0.50 to -6.00 (0.25D steps) -6.50 to -10.00 (0.50D steps)
Focus Dailies toric	Base Curve: 8.6 Diameter: 14.2 Dk/t @-3.00D 26	SPH: +0.25 to +4.00 (0.25 steps) plano to -8.00 (0.25 steps) CYL: -0.75 and -1.50 AXIS: 20, 70, 90, 110, 160, 180
Focus Dailies Progressive	Base Curve: 8.6 Diameter: 13.8 Dk/t @-3.00D 23.6	+0.25 to +5.00D (0.25D steps) -0.25 to -6.00D (0.25D steps) Incorporating a single progressive add; effective range up to +3.00



Check the website: <http://optometryportal.com/>
Follow on Twitter @optometryportal
Download the 'Optometry portal' Android App



ALCON/CIBA Contact lens parameters

Alcon/CibaVision Monthly contact lenses

Contact Lens Name	Contact Lens Specifications	Contact Lens Parameters
Air Optix Aqua	Base Curve: 8.6 Diameter: 14.2 Dk/t @-3.00D: 138	+0.25 to +6.00 (0.25D steps) -0.25 to -8.00 (0.25D steps) -8.50 to -10.00 (0.50D steps)
Air Optix Night & Day Aqua	Base Curve 8.4 & 8.6 Diameter: 13.8 Dk/t @-3.00D 175	+0.25 to +6.00 (0.25D steps) -0.25 to -8.00 (0.25D steps) -8.50 to -10.00 (0.50D steps) Plano also available
Air Optix for Astigmatism	Base Curve: 8.7 Diameter: 14.5 Dk/t (@-3.00D/-1.25x180) 108	+0.25 to +6.00 (0.25D steps) -0.25 to -6.00 (0.25D steps) -6.50 to -10.00 (0.50D steps) Cyl: -0.75, -1.25, -1.75, -2.25 (Full Circle in 10° steps)
Air Optix Aqua Multifocal	Base Curve: 8.6 Diameter: 14.2 Dk/t @-3.00D 138	+0.25 to +6.00 (0.25D steps) +0.25 to -10.00 (0.25D steps) Add Powers: LO, MED, HI



Check the website: <http://optometryportal.com/>
Follow on Twitter @optometryportal
Download the 'Optometry portal' Android App