



Sauflon Contact Lens Parameters

Sauflon Daily disposable contact lenses

Contact Lens Name	Contact Lens Specifications	Contact Lens powers	Astigmatism powers and axis available
Clariti 1 day	Base Curve: 8.6 Diameter: 14.1 Dk/t @-3.00D 86	-0.50D to -6.00D (0.25D steps) -6.50D to -10.00D (0.50D steps) +0.50D to +6.00D (0.25D steps) +6.50D to +8.00D(0.50D steps)	
Clariti 1 day toric	Base Curve: 8.6 Diameter: 14.3 Dk/t @-3.00D 57	plano to -6.00D (0.25D steps) -6.50D to -8.00D (0.50D steps) +0.25D to +4.00D (0.25D steps)	CYL: -0.75, -1.25 AXIS: 20, 70, 90, 110, 160, 180
Clariti 1 day multifocal	Base Curve: 8.6 Diameter: 14.1 Dk/t @-3.00D 86	-0.25D to -6.00D (0.25D steps) +0.25D to +5.00D (0.25D steps) LOW: up to +2.25D HIGH: +2.25D to +3.00D	



Check the website: <http://optometryportal.com/>
Follow on Twitter @optometryportal
Download the 'Optometry portal' Android App



Sauflon Contact Lens Parameters

Sauflon Monthly contact lenses

Contact Lens Name	Contact Lens Specifications	Contact Lens powers	Astigmatism powers and axis available
Clariti (monthly)	Base Curve: 8.4 & 8.8 Diameter: 14.1 Dk/t @-3.00D 86	-0.25 to -8.00 (0.25DS steps) -8.50 to -10.00 (0.50DS steps) +0.25 to +6.00 (0.25DS steps) +6.50 to +8.00 (0.50DS steps)	
Clariti Toric (monthly)	Base Curve: 8.7 Diameter: 14.4 Dk/t @-3.00D 57	Plano to -8.00 (0.25DS steps) -8.50 to -10.00 (0.50DS steps) +0.25 to +6.00 (0.25DS steps) +6.50 to +10.00(0.50DS steps)	CYL: -0.75, -1.25, -1.75, -2.25 AXIS: (Full Circle in 10°steps)
Clariti XR Toric (monthly)	Base Curve: 8.7 Diameter: 14.4 Dk/t @-3.00D 57	Plano to -8.00 (0.25DS steps) -8.50 to -10.00 (0.50DS steps) +0.25 to +6.00 (0.25DS steps) +6.50 to +10.00 (0.50DS steps)	CYL: -2.75 to -5.75 AXIS: (Full Circle in 5°steps)
Clariti Multifocal (monthly)	Base Curve: 8.7 Diameter: 14.2 Dk/t @-3.00D 86	+6.00 to -6.00 (0.25DS steps) -6.50 to -8.00 (0.50DS steps) LOW: up to +2.25D HIGH: +2.25D to +3.00D	



Check the website: <http://optometryportal.com/>
 Follow on Twitter @optometryportal
 Download the 'Optometry portal' Android App