



Johnson & Johnson Contact lens parameters

Johnson & Johnson Daily Disposable contact lenses

Contact Lens Name	Contact Lens Specifications	Contact Lens Parameters
1-day Acuvue Moist	Base Curve: 8.5&9.0 Diameter: 14.2 Dk/t @-3.00D 25.5	-0.50D to -6.00D (0.25D steps) -6.50D to -12.00D (0.50D steps) +0.50D to +6.00D (0.25D steps)
1-day Acuvue Moist for Astigmatism	Base Curve: 8.5 Diameter: 14.5 mm Dk/t @-3.00D 23.7	Plano to -6.00D (0.25D steps) Cyl powers: -0.75, -1.25, -1.75 AXIS: 10°, 20°, 60°, 70°, 80°, 90°, 100°, 110°, 120°, 160°, 170°, 180° For -2.25 Cyl limited axis available at 20°, 90°, 160°, 180° High minus power range: -6.50D to -9.00D (0.50D steps) CYL: -0.75, -1.25, -1.75 AXIS: 10°, 20°, 60°, 70°, 80°, 90°, 100°, 110°, 120°, 160°, 170°, 180° For -2.25 Cyl limited axis available at 20°, 90°, 160°, 180° Plus power range: +0.25D to +4.00D (0.25D steps) CYL: -0.75, -1.25, -1.75 Axis: 20°, 70°, 90°, 110°, 160°, 180°
1-day Acuvue Trueye	Base Curve: 8.5&9.0 Diameter: 14.2 Dk/t @-3.00D 118	-0.50D to -6.00D (0.25D steps) -6.50D to -12.00D (0.50D steps) +0.50D to +6.00D (0.25D steps)



Check the website:
<http://optometryportal.com/>
Follow on Twitter @optometryportal
Download the 'Optometry portal' Android App



Johnson & Johnson Contact lens parameters

Johnson & Johnson Monthly contact lenses 1/2

Contact Lens Name	Contact Lens Specifications	Minus Contact Lens Powers	Plus contact Lens powers	Astigmatism powers and axis available
Acuvue Advance	Base Curve: 8.3 & 8.7 Diameter: 14.0 Dk/t @-3.00D 85.7	-0.50 to -6.00 (0.25D steps) -6.50 to -12.00 (0.50D steps)	+0.50 to +6.00 (0.25D steps) +6.50 to +8.00 (0.50D steps)	
Acuvue Advance for Astigmatism	Base Curve: 8.6 Diameter: 14.5 Dk/t @-3.00D 85.7	Low minus power range: Plano to -6.00D (0.25D steps) High minus power range: -6.50D to -9.00D (0.50D steps)	Plus power range: +0.25D to +6.00D (0.25D steps)	Low minus cyl range: CYL: -0.75, -1.25, -1.75 Full Circle (10° steps) For -2.25 Cyl Axis available: 10°, 20°, 70°, 80°, 90°, 100°, 110°, 160°, 170°, 180° High minus cyl range: CYL: -1.25, -1.75 AXIS: 10°, 20°, 70°, 80°, 90°, 100°, 110°, 160°, 170°, 180° For -2.25 Cyl Axis available: 10°, 20°, 70°, 80°, 90°, 100°, 110°, 160°, 170°, 180° Plus power cyl range: CYL: -0.75, -1.25, -1.75 AXIS: 10°, 20°, 70°, 80°, 90°, 100°, 110°, 160°, 170°, 180°
Acuvue Oasys	Base Curve: 8.4 & 8.8, Diameter: 14.0 Dk/t @-3.00D 147	-0.50D to -6.00D (0.25D steps) -6.50D to -12.00D (0.50D steps)	+0.50D to +6.00D (0.25D steps) +6.50D to +8.00D (0.50D steps)	



Check the website:
<http://optometryportal.com/>
 Follow on Twitter @optometryportal
 Download the 'Optometry portal' Android App



Johnson & Johnson Contact lens parameters

Johnson & Johnson Monthly contact lenses 2/2

Contact Lens Name	Contact Lens Specifications	Minus Contact Lens Powers	Plus contact Lens powers	Astigmatism powers and axis available
Acuvue Oasys for Astigmatism	Base Curve: 8.6 Diameter: 14.5 Dk/t @-3.00D 129	Low minus power range: Plano to -6.00D (0.25D steps) High minus power range: -6.50D to -9.00D (0.50D steps)	Plus power range: +0.25D to +6.00D (0.25D steps)	Low minus cyl range: CYL: -0.75, -1.25, -1.75 Full Circle (10° steps) For -2.25 Cyl Axis available: 10°, 20°, 70°, 80°, 90°, 100°, 110°, 160°, 170°, 180° High minus cyl range: CYL: -0.75, -1.25, -1.75 Full Circle (10° steps) For -2.25 Cyl Axis available: 10°, 20°, 90°, 160°, 170°, 180° Plus power cyl range: CYL: -0.75, -1.25, -1.75 AXIS: 10°, 20°, 70°, 80°, 90°, 100°, 110°, 160°, 170°, 180° For -2.25 Cyl Axis available: 10°, 20°, 90°, 160°, 170°, 180°
Acuvue Oasys for Presbyopia	Base Curve: 8.4 Diameter: 14.3 Dk/t @-3.00D 147	plano to -9.00D (0.25D steps) ADDs: Low (+0.75D to +1.25D), Mid (+1.50D to +1.75D) High (+2.00D to +2.50D)	+0.25 to +6.00D (0.25D steps) ADDs: Low (+0.75D to +1.25D) Mid (+1.50D to +1.75D) High (+2.00D to +2.50D)	



Check the website:
<http://optometryportal.com/>
 Follow on Twitter @optometryportal
 Download the 'Optometry portal' Android App